

Product Specification Summary

SCA2

Catalog Number 40-2038-11
Product Name SCA2 Genemer Kit; 1 Kit
Size 1 Kit

Genemer™ genotyping kits are optimized components for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Infomation <http://www.genelink.com/geneprodsite/product.asp?p=801>



Product Manual http://www.genelink.com/Literature/ps/M40-2038-11_V4.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|--------|
| SCA2 22 repeat (19 CAG + 3 CAA) GScan & Genemer Control DNA | 40-2038-01 | 500 ng |
| SCA2 GeneProber unlabeled probe; 500 ng | 40-2038-40 | 500 ng |
| SCA2 Genemer Kit; 1 Kit | 40-2038-11 | 1 Kit |