



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-5020-25
Product Name	Hybwash A, Hyb Wash Solution Conc. (20X SSC); 250 mL
Size	250 mL
Description	Chemiluminescent Detection
Component/Note	Concentrate
Component/Note	Hybridization Wash Solution

Hybwash A, Hybridization Wash Solution Concentrate (20X SSC)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=83>



Product Manual http://www.genelink.com/Literature/ps/Dig_Southern_Reagents_V3.4.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL