



# Product Specification Summary

## GeneProber™ Unlabeled Probes

Catalog Number 40-2026-38  
Product Name Myotonic Dystrophy GeneProber GLDM3 unlabeled Probe; 500 ng  
Size 500 ng

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=870>



Product Manual [http://www.genelink.com/Literature/ps/M40-2026-40\\_Ver5.2.pdf](http://www.genelink.com/Literature/ps/M40-2026-40_Ver5.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
SCA2 GeneProber unlabeled probe; 500 ng	40-2038-40	500 ng
SCA3 GeneProber unlabeled probe; 500 ng	40-2039-40	500 ng