



Product Specification Summary

Genetic Tools & Reagents

| | |
|----------------|--|
| Catalog Number | 08-0021-02 |
| Product Name | Omni-Array™ mRN Antisense strand amplification kit |
| Size | 1 Kit |
| Description | Omni-Array mRNA Antisense strand amplification kit |

The Omni-Array system amplifies mRNA from as little as 2 ng. The kit offers the user two protocols for amplification of antisense strand RNA depending on the initial amount of total RNA present. A single round protocol is sufficient to generate 10 µg of antisense strand RNA from 100 ng of total RNA. When the initial amount of total RNA is less than 5 ng, a two round amplification protocol is recommended.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=898>



Product Manual <http://www.genelink.com/Literature/ps/M08-0021-02 Ver4.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit |
| Omni-Array™ Sense Strand mRNA Amplification Kit | 08-0011-02 | 1 Kit |
| Omni-Marker™ Universal unlabeled; 100 µL | 40-3005-01 | 100 µL |