



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number 40-3061-01  
Product Name GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL  
Size 100 uL

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=906>



Product Manual <http://www.genelink.com/Literature/ps/PS40-3062-XX.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product  | Catalog No | Size   |
|--|------------|--------|
| Omni-Array™ mRN Antisense strand amplification kit | 08-0021-02 | 1 Kit  |
| Omni-Array™ Sense Strand mRNA Amplification Kit    | 08-0011-02 | 1 Kit  |
| Omni-Marker™ Universal unlabeled; 100 uL           | 40-3005-01 | 100 uL |