



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5023-20
Product Name	Lumisol II, Hybridization Solution; 200 mL
Size	200 mL
Component/Note	Ready-to-use hybridization solutions specifically formulated for non-radioactive chemiluminescent hybridization.

Lumisol™ II hybridization solution can be used for all stringent types of nucleic acid blot hybridization conditions, particularly using non-radioactive digoxigenin labeled probes. Lumisol™ II hybridization solution can be used for prehybridization and hybridization. The hybridization temperature should be calculated for proper results. Duration of hybridization can be reduced to 6 hours and overnight hybridization can be used for high sensitivity requirements and for convenience.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=94>



Product Manual [http://www.genelink.com/Literature/ps/PS40-5023\\_Lumisol\\_Ver\\_1.5.pdf](http://www.genelink.com/Literature/ps/PS40-5023_Lumisol_Ver_1.5.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL