



## Product Specification Summary

### Omni-Pure™ Plasmid DNA Purification Systems

Catalog Number	40-4020-01
Product Name	Omni-Pure™ Plasmid Purification System; 100 purifications
Size	100 Purifications
Description	Omni-Pure™ Plasmid DNA Purification System

Facile and rapid purification of DNA from varied sources can be performed using the Omni-Pure™ series of DNA, RNA and plasmid purification systems. Each purification system has been formulated, optimized and designed to yield the highest purity available with the starting sample volume specially geared towards the desired downstream application. Routine mini-preps of plasmid extraction are made even easier with consistent performance. Purification can be performed with a maximum of 3 ml of cells yielding up to 20 µg of purified DNA. The convenient spin column method can be scaled up by using multiple columns and processed in less than 30 minutes. The purified DNA is of high quality suitable for all molecular biology applications including direct use in fluorescent automated sequencing methods.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=988>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-4020-01.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

### Product

### Catalog No

### Size

Omni-Pure™ Plasmid Purification System; 100 purifications

40-4020-01

100 Purifications

Omni-Pure™ Plasmid Purification System; 500 purifications

40-4020-05

500 Purifications