



[Directions to Gene Link](#)
[190 Saw Mill River Road, Hawthorne, NY 10532](#)
[\(914\) 769-1192](#)

From Long Island and Queens

From Whitestone Bridge or Throgs Neck Bridge, continue straight onto the Hutchinson River Parkway to exit 13-Cross County Parkway West. At end of Cross County Parkway, bear right onto Saw Mill Parkway North. Proceed to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North (Saw Mill River Road). Proceed about 1 mile and make a right turn on Belmont Road. Driveway ~50 ft to the left.

From George Washington Bridge

Take I-87 North-Major Deegan Expressway through Yonkers to Saw Mill Parkway North. Proceed to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North (Saw Mill River Road). Proceed about 1 mile and make a right turn on Belmont Road. Driveway ~50 ft to the left.

From Yonkers and Upper Westchester

Saw Mill Parkway to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North (Saw Mill River Road). Proceed about 1 mile and make a right turn on Belmont Road. Driveway ~50 ft to the left.

From Connecticut

Follow Cross Westchester Expressway/287 West to exit 2-9A North. Follow for 4 miles. Make a right turn on Belmont Road. Driveway ~50 ft to the left.

From West Side Manhattan

From West Side take West Side Highway which becomes Henry Hudson Parkway which becomes Saw Mill Parkway to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North (Saw Mill River Road). Proceed about 1 mile and make a right turn on Belmont Road. Driveway ~50 ft to the left.

From Upstate New York/Tappan Zee Bridge

Off Bridge take exit 8-Cross Westchester Expressway/287 East to exit 1-Saw Mill Parkway North to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North (Saw Mill River Road). Proceed about 1 mile and make a right turn on Belmont Road. Driveway ~50 ft to the left.

From East Side Manhattan

From East Side take FDR Drive to Triboro Bridge to Bruckner Expressway to Bronx River Parkway North. Continue to Taconic/Sprain Brook Parkways North. Take route 141/Hawthorne exit. Turn left at end of ramp and left again at stop sign onto route 141. Follow signs to route 9A South. From 9A South, 190 Saw Mill River Road, Gene Link building on left side. Driveway from Belmont Road across from Green valley Nursery.

From Upstate New York/East of Hudson River

Take Taconic State Parkway South to route 100 exit/Briarcliff Manor. Turn left at traffic light onto 100 South. Follow to 9A South (Saw Mill River Road) for about 2 miles. 190 Saw Mill River Road, Gene Link building on left side. Driveway from Belmont Road across from Green valley Nursery.

