



Product Specifications


Custom Oligo Synthesis, antisense oligos, RNA oligos, chimeric oligos, Fluorescent dyes, Affinity Ligands, Spacers & Linkers, Duplex Stabilizers, Minor bases, labeled oligos, Molecular Beacons, siRNA, phosphonates Locked Nucleic Acids (LNA); 2'-5' linked Oligos

Oligo Modifications

For research use only. Not for use in diagnostic procedures for clinical purposes.

Spacer C6

Category	Spacers
Modification Code	SpC6
Reference Catalog Number	26-6945
5 Prime	Y
3 Prime	Y
Internal	Y
Molecular Weight(mw)	180.14



Spacer C6 phosphoramidite is used to incorporate a six carbon spacer arm into oligonucleotides. Spacer phosphoramidites can be added sequentially to create longer spacer arms. This C6 spacer amidite has an aliphatic linker.