



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	08-0021-02
Product Name	Omni-Array™ mRN Antisense strand amplification kit
Size	1 Kit
Description	Omni-Array mRNA Antisense strand amplification kit

The Omni-Array system amplifies mRNA from as little as 2 ng. The kit offers the user two protocols for amplification of antisense strand RNA depending on the initial amount of total RNA present. A single round protocol is sufficient to generate 10 µg of antisense strand RNA from 100 ng of total RNA. When the initial amount of total RNA is less than 5 ng, a two round amplification protocol is recommended.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=898>



Product Manual <http://www.genelink.com/Literature/ps/M08-0021-02 Ver4.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 µL	40-3005-01	100 µL