



## Product Specification Summary

### Sequencing Primers

Catalog Number 26-3000-04  
Product Name M13/pUC Reverse (-48) 24mer; 25 ug  
Size 25 ug

All M13/pUC primers are complementary to sequence in the M13 and pUC vectors (1). These primers will yield sequencing data for inserts cloned in these vectors. In addition to these vectors sequencing data will be obtained for all templates/vectors containing the lacZ sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis. 1. Messing J. (1983) Meth. Enzymol. 101:20.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=99>



Product Manual [http://www.genelink.com/Literature/ps/M13seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/M13seq_primer_PSV3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	25 ug