



## Product Specification Summary

### Sequencing Primers

Catalog Number 26-3000-07  
Product Name SP6 Promoter primer 24mer; 25 ug  
Size 25 ug

This SP6 RNA polymerase primer is complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III T7 RNA polymerase promoter sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=152>



Product Manual [http://www.genelink.com/Literature/ps/T7\\_T3\\_seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	25 ug