



## Product Specification Summary

### Sequencing Primers

Catalog Number 26-3000-08  
Product Name Lambda gt11 forward primer 24mer; 25 ug  
Size 25 ug

This lambda gt11 forward primer is upstream of the unique Eco RI cloning site present in the phage imm434 gene region. This primer can be used to sequence inserts, and can be used with the Lambda gt10 reverse primer (catalog # 26-3000-11) to amplify the insert using polymerase chain reaction. The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=90>



Product Manual [http://www.genelink.com/Literature/ps/Lambda\\_seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/Lambda_seq_primer_PSV3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

| Product            | Catalog No | Size  |
|--------------------|------------|-------|
| M13 Rev Seq; 25 ug | 26-3000-52 | 25 ug |
| M13 For Seq; 25 ug | 26-3000-51 | 25 ug |
| T7 Seq Long        | 26-3000-55 | 25 ug |