



Product Specification Summary

Sequencing Primers

| | |
|----------------|---|
| Catalog Number | 26-3000-13 |
| Product Name | T7 Universal Primer (20 mer) |
| Size | 25 ug |
| Description | T7 Universal Primer (20 mer) |
| Component/Note | Sequence 5'-TAATACGACTCACTATAGGG-3' |
| Component/Note | supplied as a lyophilized powder, after reconstitution store at -20oC. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis. Same sequence with Catalog#: 26-3000-55 |

The T7, T3 and SP6 RNA polymerase primers are complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III RNA polymerase promoter sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20oC. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17394>

Product Manual  http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf

Product MSDS  <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|-------|
| M13 Rev Seq; 25 ug | 26-3000-52 | 25 ug |
| M13 For Seq; 25 ug | 26-3000-51 | 25 ug |
| T7 Seq Long | 26-3000-55 | 25 ug |
| T7 Universal Primer (20 mer) | 26-3000-13 | 25 ug |
| T7 terminator | 26-3000-56 | 25 ug |
| Lambda gt10 forward primer 21mer; 25 ug | 26-3000-10 | 25 ug |
| Lambda gt10 reverse primer 24mer; 25 ug | 26-3000-11 | 25 ug |
| Lambda gt11 forward primer 24mer; 25 ug | 26-3000-08 | 25 ug |
| Lambda gt11 reverse primer 24mer; 25 ug | 26-3000-09 | 25 ug |
| M13/pUC (-20) 17mer; 25 ug | 26-3000-01 | 25 ug |
| M13/pUC (-40) 17mer; 25 ug | 26-3000-03 | 25 ug |
| M13/pUC Reverse (-24) 16mer; 25 ug | 26-3000-02 | 25 ug |
| M13/pUC Reverse (-48) 24mer; 25 ug | 26-3000-04 | 25 ug |
| SP6 Promoter primer 24mer; 25 ug | 26-3000-07 | 25 ug |
| SP6 Universal 18mer; 25 ug | 26-3000-12 | 25 ug |