



Product Specification Summary

Sequencing Primers

| | |
|----------------|---------------------------------------|
| Catalog Number | 26-3000-21 |
| Product Name | T7- Bluescript promoter primer; 25 ug |
| Size | 25 ug |
| Description | 22 mer |
| Component/Note | GTAATACGACTCACTATAGGGC |
| Storage | Store at -20°C |

This T7 RNA polymerase primer is complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III T7 RNA polymerase promoter sequences. This primer in addition contains sequence complementary to Stratagenes's pBluescript plasmid flanking the T7 sequence. The product is supplied as a lyophilized powder, after reconstitution store at -20oC. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=161>



Product Manual

http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|-------|
| M13 Rev Seq; 25 ug | 26-3000-52 | 25 ug |
| M13 For Seq; 25 ug | 26-3000-51 | 25 ug |
| T7 Seq Long | 26-3000-55 | 25 ug |
| T7 Universal Primer (20 mer) | 26-3000-13 | 25 ug |
| T7 terminator | 26-3000-56 | 25 ug |
| Lambda gt10 forward primer 21mer; 25 ug | 26-3000-10 | 25 ug |
| Lambda gt10 reverse primer 24mer; 25 ug | 26-3000-11 | 25 ug |
| Lambda gt11 forward primer 24mer; 25 ug | 26-3000-08 | 25 ug |
| Lambda gt11 reverse primer 24mer; 25 ug | 26-3000-09 | 25 ug |
| M13/pUC (-20) 17mer; 25 ug | 26-3000-01 | 25 ug |
| M13/pUC (-40) 17mer; 25 ug | 26-3000-03 | 25 ug |
| M13/pUC Reverse (-24) 16mer; 25 ug | 26-3000-02 | 25 ug |
| M13/pUC Reverse (-48) 24mer; 25 ug | 26-3000-04 | 25 ug |
| SP6 Promoter primer 24mer; 25 ug | 26-3000-07 | 25 ug |
| SP6 Universal 18mer; 25 ug | 26-3000-12 | 25 ug |