



## Product Specification Summary

### Sequencing Primers

Catalog Number	26-3000-22
Product Name	T3- Bluescript promoter primer; 25 ug
Size	25 ug
Description	20 mer
Component/Note	AATTAACCCTCACTAAAGGG
Storage	Store at -20°C

This T3 RNA polymerase primer is complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III T3 RNA polymerase promoter sequences. This primer in addition contains sequence complementary to Stratagenes's pBluescript plasmid flanking the T3 sequence. The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=159>



Product Manual

[http://www.genelink.com/Literature/ps/T7\\_T3\\_seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/T7_T3_seq_primer_PSV3.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	25 ug
T7 Universal Primer (20 mer)	26-3000-13	25 ug
T7 terminator	26-3000-56	25 ug
Lambda gt10 forward primer 21mer; 25 ug	26-3000-10	25 ug
Lambda gt10 reverse primer 24mer; 25 ug	26-3000-11	25 ug
Lambda gt11 forward primer 24mer; 25 ug	26-3000-08	25 ug
Lambda gt11 reverse primer 24mer; 25 ug	26-3000-09	25 ug
M13/pUC (-20) 17mer; 25 ug	26-3000-01	25 ug
M13/pUC (-40) 17mer; 25 ug	26-3000-03	25 ug
M13/pUC Reverse (-24) 16mer; 25 ug	26-3000-02	25 ug
M13/pUC Reverse (-48) 24mer; 25 ug	26-3000-04	25 ug
SP6 Promoter primer 24mer; 25 ug	26-3000-07	25 ug
SP6 Universal 18mer; 25 ug	26-3000-12	25 ug