



## Product Specification Summary

### Sequencing Primers

Catalog Number	26-3000-52
Product Name	M13 Rev Seq; 25 ug
Size	25 ug
Description	M13 Sequencing Primer
Component/Note	4.83 nmols
Component/Note	5'-CAGGAAACAGCTATGAC-3'

All M13/pUC primers are complementary to sequence in the M13 and pUC vectors (1). These primers will yield sequencing data for inserts cloned in these vectors. In addition to these vectors sequencing data will be obtained for all templates/vectors containing the lacZ sequences. The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis. 1. Messing J. (1983) Meth. Enzymol. 101:20.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17145>



Product Manual

[http://www.genelink.com/Literature/ps/M13seq\\_primer\\_PSV3.1.pdf](http://www.genelink.com/Literature/ps/M13seq_primer_PSV3.1.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
M13 Rev Seq; 25 ug	26-3000-52	25 ug
M13 For Seq; 25 ug	26-3000-51	25 ug
T7 Seq Long	26-3000-55	25 ug
T7 Universal Primer (20 mer)	26-3000-13	25 ug
T7 terminator	26-3000-56	25 ug
Lambda gt10 forward primer 21mer; 25 ug	26-3000-10	25 ug
Lambda gt10 reverse primer 24mer; 25 ug	26-3000-11	25 ug
Lambda gt11 forward primer 24mer; 25 ug	26-3000-08	25 ug
Lambda gt11 reverse primer 24mer; 25 ug	26-3000-09	25 ug
M13/pUC (-20) 17mer; 25 ug	26-3000-01	25 ug
M13/pUC (-40) 17mer; 25 ug	26-3000-03	25 ug
M13/pUC Reverse (-24) 16mer; 25 ug	26-3000-02	25 ug
M13/pUC Reverse (-48) 24mer; 25 ug	26-3000-04	25 ug
SP6 Promoter primer 24mer; 25 ug	26-3000-07	25 ug
SP6 Universal 18mer; 25 ug	26-3000-12	25 ug