



Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=61>



Product Manual

<http://www.genelink.com/Literature/ps/ps26-3100-03.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Eco RI/Bam HI Adaptor; 20 ug	26-3100-22	20 ug
Sal I/Bam HI Adaptor; 20 ug	26-3100-23	20 ug
Hind III/Eco RI Adaptor; 20 ug	26-3100-20	20 ug
Bam HI/Bgl II Adaptor; 20 ug	26-3100-21	20 ug
XhoI/Hind III Adaptor	26-3100-24	20 ug
Bgl II /Not I Adaptor; 20 ug	26-3100-19	20 ug
Eco RI/ Not I Adaptor; 20 ug	26-3100-02	20 ug
Eco RI/ Sma I Adaptor; 20 ug	26-3100-03	20 ug
Not I/ Kpn I adaptor; 20 ug	26-3100-14	20 ug
Sac I/Not I Adaptor; 20 ug	26-3100-10	20 ug
Sph I/Not I adaptor; 20 ug	26-3100-13	20 ug
Sph/Spe I Adaptor; 20 ug	26-3100-09	20 ug
Xba I/EcoR I adaptor; 20 ug	26-3100-15	20 ug
Xho I/ Eco RI Adaptor; 20 ug	26-3100-05	20 ug
Xho I/ Not I Adaptor; 20 ug	26-3100-06	20 ug