

# Product Specification Summary

Catalog Number	26-4000-10
Product Name	Random Hexamer Phosphorylated pd(N)6; 50 ug
Size	50 ug
Description	5'pNNNNNN
Component/Note	mw: 1872
Component/Note	nmol: ~26
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=138>



Product Manual <http://www.genelink.com/Literature/ps/PS26-4000-10V1.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
5'-Dig Random Hexamer; 25 ug	26-4000-81	25 ug
5'-Dig Random Nonamer; 25 ug	26-4000-84	25 ug

5'-Dig Random Octamer; 25 ug  
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS26-4000-10.pdf Print Date Version : October 20, 2018 Page 1 of 1