

Product Specification Summary

Fluorescent Molecular Probes

Catalog Number 26-8122-10
Product Name Fluorophore-Quencher Probe, 5'-Hex/3'-Quencher; 600 nmols yield
Size 600 nmols
Description Fluorophore-Quencher Probe, 5'-Hex/3'-Quencher; 600 nmols yield
Scan the QR Code or visit the following links
Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17419>

Product Manual http://www.genelink.com/Literature/ps/Fluorescent_oligo_recons_V1.pdf

Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>

Related Products

Product	Catalog No	Size
SmartBase™ Fluorescent Molecular Probes	SB-FMP	1 each
QPCR_Manual	M_QPCR	1 each
5'-Fam/3'-Quencher 10 umol scale; 600 nmol yield	26-8121-10	600 nmols