



# Product Specification Summary

## RNA Interference

Catalog Number	27-6402-01
Product Name	Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit
Size	1 Kit
Description	Guaranteed RNAi Explorer kit with Fluorescein/Dabcyl Molecular Beacon

3+1 siRNAs and a detection probe. You simply supply the accession number or sequence for the gene to down regulate. Gene Link designs siRNAs using the latest published criteria and proprietary algorithms. This kit contains three target siRNAs, a negative control and a detection probe of your choice OR 6 custom oligos\*.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=963>



Product Manual [http://www.genelink.com/Literature/ps/RNAi\\_Guar\\_PS\\_Ver\\_4.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver_4.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit