



Product Specification Summary

RNA Interference

Catalog Number 27-6402-06
Product Name Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit
Size 1 Kit
Description Guaranteed RNAi Explorer kit with 6 oligos up to 30mer at 50 nmol scale

3+1 siRNAs and a detection probe. You simply supply the accession number or sequence for the gene to down regulate. Gene Link designs siRNAs using the latest published criteria and proprietary algorithms. This kit contains three target siRNAs, a negative control and a detection probe of your choice OR 6 custom oligos*.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=965>



Product Manual [http://www.genelink.com/Literature/ps/RNAi Guar PS Ver 4.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver_4.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each | 27-6402-02 | 1 Kit |
| Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit | 27-6402-06 | 1 Kit |