



## Product Specification Summary

### RNA Interference

Catalog Number	27-6416-05
Product Name	Human GAPDH control siRNA duplex unlabeled; 5 nmols
Size	5 nmols
Description	Positive Control
Component/Note	Unlabeled
Component/Note	GAPDH siRNA Duplex
Storage	Store at -20°C

A negative control should also be included when performing siRNA experiments to eliminate the possibility of nonspecific silencing effects. For this purpose, Gene Link provides a siRNA sequence that has been shown to have no effect on gene silencing. The negative control siRNA is also available labeled with fluorescein, which allows transfection efficiency to be tracked by fluorescent microscopy.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=11659>



Product Manual

[http://www.genelink.com/Literature/ps/RNAi\\_Controls\\_COA\\_PS\\_V5.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Controls_COA_PS_V5.1.pdf)



Product MSDS

[http://www.genelink.com/Literature/ps/RNAi\\_Controls\\_PS\\_Ver\\_4.3.pdf](http://www.genelink.com/Literature/ps/RNAi_Controls_PS_Ver_4.3.pdf)



## Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit
Negative Control non-silencing FL labeled; 20 nmols	27-6411-20FL	20 nmols
Human Vimentin custom siRNA duplex; 20 nmols	27-6412-20	20 nmols
Mouse Beta Tubulin control siRNA duplex; 20 nmols	27-6414-20	20 nmols
Rat Chromogranin-A control siRNA duplex; 20 nmols	27-6415-20	20 nmols
Human Vimentin control siRNA duplex FL label; 10 nmols	27-6412-10FL	10 nmols
Human Beta Actin control siRNA duplex FL; 10 nmols	27-6413-10FL	10 nmols
Rat Chromogranin-A control siRNA duplex; 10 nmols	27-6415-10FL	10 nmols
Custom siRNA duplex standard purity; 2 nmols	27-6436-20	2 nmols
Custom siRNA duplex standard purity; 10 nmols	27-6436-10	10 nmols
Custom siRNA duplex standard purity; 50 nmols	27-6436-05	50 nmols
Custom siRNA duplex standard purity; 100 nmols	27-6436-01	100 nmols
Custom siRNA duplex Reverse Phase Purified; 2 nmols	27-6437-20	2 nmols