



Product Specification Summary

RNA Interference

Catalog Number	27-6430-05
Product Name	Pre-designed 3 siRNA standard purity; 5 nmols
Size	5 nmols
Description	Pre-designed siRNA (3+) ; standard purity, 5 nmol
Component/Note	Pre-designed siRNA (3+) ; standard purity, 5 nmol
Component/Note	Pre-designed siRNA (3+) ; standard purity, 5 nmol
Pre-designed siRNA (3+) ; standard purity, 5 nmol	

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=16982>



Product Manual http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver5.2.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit