



# Product Specification Summary

## RNA Interference

Catalog Number 27-6430-05  
Product Name Pre-designed 3 siRNA standard purity; 5 nmols  
Size 5 nmols  
Description Pre-designed siRNA (3+) ; standard purity, 5 nmol  
Component/Note Pre-designed siRNA (3+) ; standard purity, 5 nmol  
Component/Note Pre-designed siRNA (3+) ; standard purity, 5 nmol  
Pre-designed siRNA (3+) ; standard purity, 5 nmol

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=16982>



Product Manual [http://www.genelink.com/Literature/ps/RNAi\\_Guar\\_PS\\_Ver5.2.pdf](http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver5.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit