



Product Specification Summary

RNA Interference

| | |
|---|---|
| Catalog Number | 27-6430-05 |
| Product Name | Pre-designed 3 siRNA standard purity; 5 nmols |
| Size | 5 nmols |
| Description | Pre-designed siRNA (3+) ; standard purity, 5 nmol |
| Component/Note | Pre-designed siRNA (3+) ; standard purity, 5 nmol |
| Component/Note | Pre-designed siRNA (3+) ; standard purity, 5 nmol |
| Pre-designed siRNA (3+) ; standard purity, 5 nmol | |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=16982>



Product Manual http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver5.2.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each | 27-6402-02 | 1 Kit |
| Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit | 27-6402-06 | 1 Kit |