



Product Specification Summary

RNA Interference

Catalog Number	27-6437-20
Product Name	Custom siRNA duplex Reverse Phase Purified; 2 nmols
Size	2 nmols
Description	Custom siRNA duplex Reverse Phase Purified, 2 nmols

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and μg . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at oligos@genelink.com

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17108>



Product Manual http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS27-6437-20.pdf Print Date Version : June 18, 2018

27-6402-06
Page 1 of 1

