



# Product Specification Summary

## RNA Interference

|                |   |
|----------------|---|
| Catalog Number | 27-6437-50  |
| Product Name   | Custom siRNA duplex RPC Purified; 5 nmols           |
| Size           | 5 nmols   |
| Description    | Custom siRNA duplex Reverse Phase Purified; 5 nmols |

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and  $\mu\text{g}$ . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at [oligos@genelink.com](mailto:oligos@genelink.com)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17324>



Product Manual [http://www.genelink.com/Literature/ps/siRNA\\_Custom\\_Duplex\\_PS\\_V1.1.pdf](http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product  | Catalog No | Size  |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each         | 27-6402-02 | 1 Kit |

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit  
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS27-6437-50.pdf Print Date Version : July 26, 2017

27-6402-06  
Page 1 of 1

