



# Product Specification Summary

## RNA Interference

Catalog Number 27-6438-10  
Product Name Custom siRNA duplex PAGE Purified; 10 nmols  
Size 10 nmols  
Description Custom siRNA duplex PAGE Purified, 10 nmols

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and  $\mu\text{g}$ . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at [oligos@genelink.com](mailto:oligos@genelink.com)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17113>



Product Manual [http://www.genelink.com/Literature/ps/siRNA\\_Custom\\_Duplex\\_PS\\_V1.1.pdf](http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product  | Catalog No | Size  |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each         | 27-6402-02 | 1 Kit |

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit 27-6402-06  
190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS27-6438-10.pdf Print Date Version : February 16, 2019 Page 1 of 1

