



# Product Specification Summary

## RAPD Decamer Sets

Catalog Number 40-0001-00  
Product Name GL RAPD Decamer Set 1- 30  
Size 1 set

Collection of all 30 sets of GL RAPD Decamer 01 through 30. The total number of oligos is 1200. Arbitrary sequence decamer sets aid in rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments. Each tube contains 10 nmols of lyophilized primer. This is sufficient for thousands of amplification reaction. Please see Protocol section for details. Please see the pdf file for detailed information.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=752>



Product Manual [http://www.genelink.com/Literature/ps/M40-0001-XX\\_Ver3.2.pdf](http://www.genelink.com/Literature/ps/M40-0001-XX_Ver3.2.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
GL RAPD Decamer Set 1- 30	40-0001-00	1 set
GL RAPD Decamer 10-ST Set	40-0001-10	1 set
GL RAPD Decamer 11-UV Set	40-0001-11	1 set