



## Product Specification Summary

### RAPD Decamer Sets

Catalog Number	40-0001-04
Product Name	GL RAPD Decamer 04-GH Set
Size	1 set
Description	GL RAPD Decamer 04-GH Set
Component/Note	Supplied Lyophilized
Component/Note	40 Decamers 10 nmols each
Storage	Store at -20°C

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments. Each tube contains 10 nmols of lyophilized primer. This is sufficient for thousands of amplification reaction. Please see Protocol section for details. Please see the pdf file for detailed information

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=724>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-0001-04.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
GL RAPD Decamer Set 1- 30	40-0001-00	1 set
GL RAPD Decamer 10-ST Set	40-0001-10	1 set
GL RAPD Decamer 11-UV Set	40-0001-11	1 set
GL RAPD Decamer 13-YZ Set	40-0001-13	1 set
GL RAPD Decamer 16-AEAF Set	40-0001-16	1 set
GL RAPD Decamer 17-AGAH Set	40-0001-17	1 set
GL RAPD Decamer 18-AIAJ Set	40-0001-18	1 set
GL RAPD Decamer 19-AKAL Set	40-0001-19	1 set
GL RAPD Decamer 20-AMAN Set	40-0001-20	1 set
GL RAPD Decamer 21-AOAP Set	40-0001-21	1 set
GL RAPD Decamer 22-AQAR Set	40-0001-22	1 set
GL RAPD Decamer 23-ASAT Set	40-0001-23	1 set
GL RAPD Decamer 24-AUAV Set	40-0001-24	1 set
GL RAPD Decamer 25-AWAX Set	40-0001-25	1 set
GL RAPD Decamer 26-AYAZ Set	40-0001-26	1 Set