



# Product Specification Summary

## RAPD Decamer Sets

|                |                           |
|----------------|---------------------------|
| Catalog Number | 40-0001-06                |
| Product Name   | GL RAPD Decamer 06-KL Set |
| Size           | 1 set                     |
| Description    | GL RAPD Decamer 06-KL Set |
| Component/Note | Supplied Lyophilized      |
| Component/Note | 40 Decamers 10 nmols each |

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=726>



Product Manual <http://www.genelink.com/Literature/ps/PS40-0001-06.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                   | Catalog No | Size  |
|---------------------------|------------|-------|
| GL RAPD Decamer Set 1- 30 | 40-0001-00 | 1 set |
| GL RAPD Decamer 10-ST Set | 40-0001-10 | 1 set |