



Product Specification Summary

RAPD Decamer Sets

Catalog Number	40-0001-15
Product Name	GL RAPD Decamer 15-ACAD Set
Size	1 set
Description	GL RAPD Decamer 15-ACAD Set
Component/Note	Supplied Lyophilized
Component/Note	40 Decamers 10 nmols each
Storage	Store at -20°C

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=735>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-0001-15.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
GL RAPD Decamer Set 1- 30	40-0001-00	1 set
GL RAPD Decamer 10-ST Set	40-0001-10	1 set
GL RAPD Decamer 11-UV Set	40-0001-11	1 set
GL RAPD Decamer 13-YZ Set	40-0001-13	1 set
GL RAPD Decamer 16-AEAF Set	40-0001-16	1 set
GL RAPD Decamer 17-AGAH Set	40-0001-17	1 set
GL RAPD Decamer 18-AIAJ Set	40-0001-18	1 set
GL RAPD Decamer 19-AKAL Set	40-0001-19	1 set
GL RAPD Decamer 20-AMAN Set	40-0001-20	1 set
GL RAPD Decamer 21-AOAP Set	40-0001-21	1 set
GL RAPD Decamer 22-AQAR Set	40-0001-22	1 set
GL RAPD Decamer 23-ASAT Set	40-0001-23	1 set
GL RAPD Decamer 24-AUAV Set	40-0001-24	1 set
GL RAPD Decamer 25-AWAX Set	40-0001-25	1 set
GL RAPD Decamer 26-AYAZ Set	40-0001-26	1 Set