



# Product Specification Summary

## RAPD Decamer Sets

Catalog Number 40-0001-16  
Product Name GL RAPD Decamer 16-AEAF Set  
Size 1 set  
Description GL RAPD Decamer 16-AEAF Set

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=736>



Product Manual <http://www.genelink.com/Literature/ps/PS40-0001-16.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product                   | Catalog No | Size  |
|---------------------------|------------|-------|
| GL RAPD Decamer Set 1- 30 | 40-0001-00 | 1 set |
| GL RAPD Decamer 10-ST Set | 40-0001-10 | 1 set |
| GL RAPD Decamer 11-UV Set | 40-0001-11 | 1 set |