



# Product Specification Summary

## RAPD Decamer Sets

Catalog Number 40-0001-26  
Product Name GL RAPD Decamer 26-AYAZ Set  
Size 1 Set  
Description GL RAPD Decamer 26-AYAZ Set

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=746>



Product Manual <http://www.genelink.com/Literature/ps/PS40-0001-26.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
GL RAPD Decamer Set 1- 30	40-0001-00	1 set
GL RAPD Decamer 10-ST Set	40-0001-10	1 set
GL RAPD Decamer 11-UV Set	40-0001-11	1 set