



# Product Specification Summary

## RT-PCRmer™

Catalog Number	40-1001-10
Product Name	beta actin A Human RT-PCRmer™ (F1/R1); Exon 3; 10 nmols
Size	10 nmols
Description	Human B- Actin Primer Mix
Component/Note	Hum-B-actin F1/R1
Storage	Store at -20°C

html

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=146>



Product Manual <http://www.genelink.com/Literature/ps/M40-1001-10V3.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols
beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols	40-1014-10	10 nmols