



# Product Specification Summary

## RT-PCRmer™

Catalog Number 40-1001-10  
Product Name beta actin A Human RT-PCRmer™ (F1/R1); Exon 3; 10 nmols  
Size 10 nmols  
Description Human B- Actin Primer Mix  
Component/Note Hum-B-actin F1/R1  
Storage Store at -20°C

html

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=146>



Product Manual <http://www.genelink.com/Literature/ps/M40-1001-10V3.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product   | Catalog No | Size     |
|---|------------|----------|
| beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols | 40-1023-10 | 10 nmols |
| beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols    | 40-1025-10 | 10 nmols |
| beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols     | 40-1014-10 | 10 nmols |