



## Product Specification Summary

### RT-PCRmer™

Catalog Number	40-1002-00
Product Name	beta actin control PCR mix (human & Rat); 200 uL
Size	200 uL
Description	Human, Rat and Mouse
Component/Note	Ready to use
Component/Note	Add Template &Taq Pol.
Storage	Store at -20°C

html head meta http-equiv="Content-Type" content="text/html; charset=windows-1252" meta name="GENERATOR" content="Microsoft FrontPage 4.0" meta name="ProgId" content="FrontPage.Editor.Document"> title>New Page 1/title> /head> /title> '>

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=52>



Product Manual

<http://www.genelink.com/Literature/ps/ps40-1002-00.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols
beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols	40-1014-10	10 nmols
beta actin F Human RT-PCRmer™ (F6/R2) 5' UTR- 3' ; 10 nmols	40-1024-10	10 nmols
beta actin Human TaqMan Probe (F2/R3); Exon 5; 2 nmols	40-1027-02	2 nmols
beta actin Human Molecular Beacon Probe (F2/R3B) Exon 5; 2 nmols	40-1029-02	2 nmols
Human beta actin MB RT-PCRmer™ for Molecular Beacon assay (F2/R3B); Exon 4-5	40-1028-10	10 nmol
beta actin Human RT-PCRmer™ for Taqman (F2/R3); 10 nmols	40-1026-10	10 nmols
beta actin B Human RT-PCRmer™ (F4/R3B) Exon 2-5; 10 nmols	40-1021-10	10 nmols
beta actin C Human RT-PCRmer™ (F6/R1) 5'UTR-Exon ; 10 nmols	40-1022-10	10 nmols
GAPDH B Human RT-PCRmer™ (F2/R2) RT-PCRmer; 10 nmols	40-1006-10	10 nmols
GAPDH C Human RT-PCRmer™ (F1/R2) Exon 2-7; 10 nmols	40-1007-10	10 nmols
GAPDH D Human RT-PCRmer™ (F1/R3) Exon 2-9; 10 nmols	40-1008-10	10 nmols
beta2 microglobulin Human RT-PCRmer ; 10 nmols	40-1004-10	10 nmols
Beta actin control PCR mix (human & Rat); 10 nmols	40-1002-10	10 nmols