



# Product Specification Summary

## RT-PCRmer™

Catalog Number	40-1003-10
Product Name	beta actin A Mouse RT-PCRmer™ (F1/R1) 3' Region; 10 nmols
Size	10 nmols
Description	Mouse beta actin
Component/Note	Primer Pair F1/R1
Storage	Store at -20°C

html

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=147>



Product Manual <http://www.genelink.com/Literature/ps/ps40-1003-10.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols
beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols	40-1014-10	10 nmols