



## Product Specification Summary

### RT-PCRmer™

Catalog Number	40-1006-10
Product Name	GAPDH B Human RT-PCRmer™ (F2/R2) RT-PCRmer; 10 nmols
Size	10 nmols
Description	GAPDH-B Primer Pair
Component/Note	Primer 2: GAPDH R2
Component/Note	Primer 1: GAPDH F2
Storage	Store at -20°C

RT-PCRmer™ are primer pairs for specific amplification of cDNA. Glyceraldehyde 3 phosphate dehydrogenase (G3PDH), like  $\beta$ -actin, is ubiquitously expressed and serves as a positive control for northern and other expression studies. These are generally used as controls for measuring cDNA synthesis efficiency by reverse transcription and as controls for mRNA (cDNA) quantitative expression studies. The enclosed human G3PDH RT-PCRmer™ are supplied as a lyophilized powder in aliquots of 10nmols. The 10 nmols of primer when dissolved in 500  $\mu$ l sterile water or TE will give a solution of 20  $\mu$ Molar i.e. 20 pmols/ $\mu$ l. The quantity supplied is sufficient for at least 400 regular 25  $\mu$ l PCR reaction\* for ethidium bromide stained visualization. The amplified product may also be detected by hybridization using Gene Link OligoProber™ with either a free 5' OH (Catalog No. 40-1105-02) for 32P labeling or biotinylated probe (Catalog No. 40-1105-03) for non-radioactive detection.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=951>



Product Manual

<http://www.genelink.com/Literature/ps/M40-1005-10 Ver3.1.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols
beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols	40-1014-10	10 nmols
beta actin F Human RT-PCRmer™ (F6/R2) 5' UTR- 3' ; 10 nmols	40-1024-10	10 nmols
beta actin Human TaqMan Probe (F2/R3); Exon 5; 2 nmols	40-1027-02	2 nmols
beta actin Human Molecular Beacon Probe (F2/R3B) Exon 5; 2 nmols	40-1029-02	2 nmols
Human beta actin MB RT-PCRmer™ for Molecular Beacon assay (F2/R3B); Exon 4-5	40-1028-10	10 nmol
beta actin Human RT-PCRmer™ for Taqman (F2/R3); 10 nmols	40-1026-10	10 nmols
beta actin B Human RT-PCRmer™ (F4/R3B) Exon 2-5; 10 nmols	40-1021-10	10 nmols
beta actin C Human RT-PCRmer™ (F6/R1) 5'UTR-Exon ; 10 nmols	40-1022-10	10 nmols
GAPDH B Human RT-PCRmer™ (F2/R2) RT-PCRmer; 10 nmols	40-1006-10	10 nmols
GAPDH C Human RT-PCRmer™ (F1/R2) Exon 2-7; 10 nmols	40-1007-10	10 nmols
GAPDH D Human RT-PCRmer™ (F1/R3) Exon 2-9; 10 nmols	40-1008-10	10 nmols
beta2 microglobulin Human RT-PCRmer ; 10 nmols	40-1004-10	10 nmols
Beta actin control PCR mix (human & Rat); 10 nmols	40-1002-10	10 nmols