



# Product Specification Summary

## RT-PCRmer™

Catalog Number	40-1027-02
Product Name	beta actin Human TaqMan Probe (F2/R3); Exon 5; 2 nmols
Size	2 nmols
Description	Beta actin TaqMan Probe
Component/Note	For F2/R3 Primers
Component/Note	TaqMan Probe F2/R3
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=16984>



Product Manual [http://www.genelink.com/Literature/ps/M40-1027-02\\_Ver1.1.pdf](http://www.genelink.com/Literature/ps/M40-1027-02_Ver1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols
beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols	40-1014-10	10 nmols