



# Product Specification Summary

## RT-PCRmer™

|                |  |
|----------------|--|
| Catalog Number | 40-1028-10   |
| Product Name   | Human beta actin MB RT-PCRmer™ for Molecular Beacon assay (F2/R3B); Exon 4-5 |
| Size           | 10 nmol  |
| Description    | Beta actin Primer Mix  |
| Component/Note | For MB Assay   |
| Component/Note | Hu-BA F2/R3B. Primer Mix   |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17178>



Product Manual <http://www.genelink.com/Literature/ps/M40-1001-10%20Ver2.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product   | Catalog No | Size     |
|---|------------|----------|
| beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols | 40-1023-10 | 10 nmols |
| beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols    | 40-1025-10 | 10 nmols |
| beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols     | 40-1014-10 | 10 nmols |