



Product Specification Summary

RT-PCRmer™

| | |
|----------------|--|
| Catalog Number | 40-1028-10 |
| Product Name | Human beta actin MB RT-PCRmer™ for Molecular Beacon assay (F2/R3B); Exon 4-5 |
| Size | 10 nmol |
| Description | Beta actin Primer Mix |
| Component/Note | For MB Assay |
| Component/Note | Hu-BA F2/R3B. Primer Mix |
| Storage | Store at -20°C |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17178>



Product Manual <http://www.genelink.com/Literature/ps/M40-1001-10%20Ver2.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|------------|----------|
| beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols | 40-1023-10 | 10 nmols |
| beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols | 40-1025-10 | 10 nmols |

beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols

40-1014-10