



Product Specification Summary

RT-PCRmer™

Catalog Number	40-1028-10
Product Name	Human beta actin MB RT-PCRmer™ for Molecular Beacon assay (F2/R3B); Exon 4-5
Size	10 nmol
Description	Beta actin Primer Mix
Component/Note	For MB Assay
Component/Note	Hu-BA F2/R3B. Primer Mix
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17178>



Product Manual <http://www.genelink.com/Literature/ps/M40-1001-10%20Ver2.1.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols

beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols

40-1014-10

190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS40-1028-10.pdf

Print Date Version : December 14, 2018

Page 1 of 1