



# Product Specification Summary

## RT-PCRmer™

Catalog Number	40-1029-02
Product Name	beta actin Human Molecular Beacon Probe (F2/R3B) Exon 5; 2 nmols
Size	2 nmols
Description	Beta actin MB Probe
Component/Note	For F2/R3B Primers
Component/Note	Molecular Beacon Probe
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17179>



Product Manual [http://www.genelink.com/Literature/ps/M40-1029-02\\_Ver1.1.pdf](http://www.genelink.com/Literature/ps/M40-1029-02_Ver1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols

beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7; 10 nmols

40-1014-10

190 Saw Mill River Road, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS40-1029-02.pdf Print Date Version : December 14, 2018

Page 1 of 1