



# Product Specification Summary

## Gene Detection Systems

Catalog Number	40-2001-01
Product Name	Sickle Cell Genemer HbA control DNA; 500 ng
Size	500 ng
Description	Sickle Cell HbA DNA
Component/Note	HbA Fragment
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=195>



Product Manual [http://www.genelink.com/Literature/ps/PS40-2001-0X\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/PS40-2001-0X_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit
Custom Genotyping	40-2000-01	1 Ea