



Product Specification Summary

Gene Detection Systems

| | |
|----------------|---|
| Catalog Number | 40-2001-01 |
| Product Name | Sickle Cell Genemer HbA control DNA; 500 ng |
| Size | 500 ng |
| Description | Sickle Cell HbA DNA |
| Component/Note | HbA Fragment |
| Storage | Store at -20°C |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=195>



Product Manual http://www.genelink.com/Literature/ps/PS40-2001-0X_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|-------------|-------|
| Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions | 40-2035-10K | 1 Kit |
| Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit | 40-2036-11 | 1 Kit |
| Custom Genotyping | 40-2000-01 | 1 Ea |