



Product Specification Summary

Gene Detection Systems

Catalog Number	40-2001-03
Product Name	Sickle Cell Genemer HbC control DNA; 500 ng
Size	500 ng
Description	Sickle Cell HbC DNA
Component/Note	HbC Fragment
Storage	Store at -20°C

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=197>



Product Manual http://www.genelink.com/Literature/ps/PS40-2001-0X_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit
Custom Genotyping	40-2000-01	1 Ea