



# Product Specification Summary

## Gene Detection Systems

Catalog Number	40-2004-20
Product Name	GLFX Genemer Kit for Radioactive Detection; 1 run
Size	1 run
Description	GLFX Genemer Kit for Radioactive Detection
Component/Note	GLFX Genemer Kit for Radioactive Detection [Old Number 40-3103-00 ]

GLFX Genemer Kit for Radioactive Detection [Old Number 40-3103-00 ] Kit for amplification and radioactive detection of Fragile X CGG triple repeat region amplified PCR products using 35S or 32P. 100 Reactions. For research use only. Not for use in diagnostic procedures for clinical purposes

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=78>



Product Manual [http://www.genelink.com/Literature/ps/M40-2004-20\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-2004-20_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit
Custom Genotyping	40-2000-01	1 Ea