



## Product Specification Summary

### Gene Detection Systems

Catalog Number	40-2004-32
Product Name	Fragile X Syndrome PCRProber Kit; 5 blots
Size	5 blots
Description	Gene Detection System
Component/Note	Contains Components A, B, C & D only
Component/Note	Amplification Components
Storage	Store at -20°C

Fragile X PCRProber AP labeled probe [Old Number 40-3101-01 ] Alkaline phosphatase labeled probe

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=76>



Product Manual

<http://www.genelink.com/Literature/ps/M40-2004-32 Ver3.1.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Factor V Leiden Genotyping Kit; 1 Kit 100 Reactions	40-2035-10K	1 Kit
Factor VIII (Hemophilia) Genemer™ Kit; 1 Kit	40-2036-11	1 Kit
Custom Genotyping	40-2000-01	1 Ea
ALS C9orf72 hexanucleotide expansion Genotyping GeneProber ALS-GL577	40-2057-41	110 uL
Factor VIII Genemer Pack (Hemophilia A); 10 nmols	40-2036-10	10 nmols
Factor V Leiden Genemer; 2 nmols	40-2035-10	2 nmols
X alphoid repeat Genemer; 10 nmoles	40-2021-10	10 nmols
Y alphoid repeat Genemer; 10 nmols	40-2022-10	10 nmols
Sickle Cell Genemer HbA control DNA; 500 ng	40-2001-01	500 ng
Sickle Cell Genemer HbS control DNA; 500 ng	40-2001-02	500 ng
Sickle Cell Genemer HbC control DNA; 500 ng	40-2001-03	500 ng
Fragile X CGG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL	40-2004-41	110 uL
Fragile X Syndrome GeneProber unlabeled Probe; 500 ng	40-2004-40	500 ng
Fragile X Syndrome PCRProber Kit; 5 blots	40-2004-32	5 blots
GLFX Genemer Kit for Radioactive Detection; 1 run	40-2004-20	1 run